

GB

## **tommee tippee® closer to nature® manual breast pump**

**Please read and retain this important information for future reference.**

Researched and developed by specialists and breastfeeding experts to express milk for storage and later use.

### **A. Parts guide** (see diagram A)

- A. Handle
- B. Pump body
- C. Bottle neck adaptor
- D. Bottle
- E. Valve (spare included)
- F. Cup with massage ripples
- G. Concertina
- H. Soft collar

**B. BEFORE USE ENSURE THE BREAST PUMP HAS BEEN DISASSEMBLED, WASHED AND STERILISED.**

### **C. DISASSEMBLING THE MANUAL BREAST PUMP** (see diagram C)

1. Unscrew bottle from bottle neck adaptor.
2. Unscrew bottle neck adaptor from pump body.
3. Pull off valve.
4. Unclip handle and pull to the side for cleaning. **IT IS NOT ESSENTIAL TO TAKE THE ASSEMBLED HANDLE APART.**

### **D. CLEANING AND STERILISING**

- Wash in warm soapy water and rinse in clean water or use a dishwasher. Do not clean with solvents or harsh chemicals.
- Do not allow product to come into contact with oil based foods (e.g. oil/tomato based sauces) as staining will occur.
- Suitable for all forms of sterilisation. Follow the steriliser manufacturer's instructions. If using the Steriliser Box follow instructions in section "Using the Steriliser Box".
- For sterilising the breast pump it is **NOT ESSENTIAL TO TAKE THE ASSEMBLED HANDLE APART.**

### **E. ASSEMBLING THE MANUAL BREAST PUMP** (see diagram E)

1. Ensure all parts of breast pump have been washed and sterilised.
2. Wash your hands thoroughly before handling the pump.
3. Push and click handle into place.
4. Push valve into place at the bottom of the pump body.
5. Attach the bottle neck adaptor to the bottle.

### **F. USING THE MANUAL BREAST PUMP**

1. Wash your hands before touching the breast pump or your breasts and avoid touching inside the breast pump or the storage pot.
2. Find a comfortable chair where you can relax
3. Hold the pump in whichever hand is most comfortable. Press the cup firmly against your breast so that it creates a seal all the way round. IF AIR CAN GET IN THERE WILL BE NO SUCTION.
4. Gently squeeze the handle and you will feel the suction. Release the handle and allow it to return to its resting position. Releasing the handle opens the valve so that the milk can flow into the bottle. Try to imitate the rhythm of your baby as he sucks. You won't need to squeeze hard – just as much as is comfortable. When you first start expressing milk, use the pump for only 3-5 minutes on each breast. Over time you will be able to build up to 10-15 minutes on each breast.
5. You can help stimulate your milk flow by gently stroking your breast.
6. Your milk will start to be released in drops and will gradually build up to a stream. First you will release your “fore” milk, a thin yellowish liquid. This will be followed by your “hind” milk which resembles skimmed milk. Don't worry that it looks thin and watery – both fore and hind milk are ideal for your baby and rich in nutrients.
7. The amount of milk you express and the time it takes will vary – don't worry, this is normal. Usually you will express 30-75mls (1-3oz) at a time, and it may take up to 30 minutes to drain full breasts. The first few times you express, the quantity may be less than this.
8. After use, take the pump apart and wash it to remove milk deposits. You do not need to sterilise your pump until you are ready to use it again.
9. Cover and store the expressed milk for later use.

**For hygiene reasons this product cannot be returned once opened unless faulty.**

#### **G. Helpful Advice**

- The best time to use your breast pump is after the first feed of the day (unless otherwise advised by your breastfeeding advisor)
- Be comfortable and relax, e.g. watch the TV or have your favourite music playing.
- Remember you are only trying to express milk from just behind the nipple – strong suction is not necessary.
- If water gets inside the **concertina** part during sterilisation, simply lift the handle slightly and pour out the water. If the **concertina** come away from the soft collar, just push it back into place.
- There should be no pain using your breast pump. If you feel pain, try using the pump more gently. Then examine your breasts. If there are any red or hot areas contact your breastfeeding advisor as soon as possible.
- Do not continue to use the pump once your milk has stopped flowing.
- Wash all breastfeeding equipment after use and before sterilisation.

**THE CLOSER TO NATURE® MANUAL BREAST PUMP CAN ONLY BE USED WITH TOMMEE TIPPEE® CLOSER TO NATURE® BOTTLES.**

**A SEPARATE INSTRUCTION SHEET ON HOW TO CLEAN AND STERILISE THE TOMMEE TIPPEE® CLOSER TO NATURE® BOTTLE IS ENCLOSED. PLEASE READ AND RETAIN THIS IMPORTANT INFORMATION FOR FUTURE REFERENCE.**

---

#### H. USING THE STERILISER BOX

The Sterilising box has been specifically designed to microwave steam sterilise all of your breastfeeding equipment quickly and efficiently.

1. Pour 100ml of water into sterilising box
2. Load upturned bottle, disassembled breast pump and any accessories.
3. Securely fit box lid
4. Open valve, on lid of box
5. Place steriliser box in microwave
6. Follow instructions according to microwave output
7. CAUTION: To avoid risk of burning, always ensure box is cool to touch before lifting out.
8. Immediately fit valve on lid to make box airtight and sterile.

| Wattage     | Full Power | Cooling Time |
|-------------|------------|--------------|
| Under 800 W | 6 minutes  | 3 minutes    |
| Over 850 W  | 4 minutes  | 3 minutes    |

---

#### I. USING THE BREAST MILK STORAGE SYSTEM (see diagram I)

For use with the Tommee Tippee® closer to nature® breast pumps

1. Ensure pot is clean and sterilised.
2. Place pot inside of bottle.
3. Screw bottle onto bottle neck adaptor and continue to follow assembly of breast pump instructions.
4. When the breast pump is fully assembled and in usage, the pot will hold up to 60ml / 2oz.
5. Once you have expressed 60ml stop pumping, place lid on storage pot and remove from bottle.
6. Store pot in either the fridge or freezer. See storage guide below.
7. If you want to continue expressing, either insert another clean and sterilised pot or express directly into the bottle of the breast pump.

#### J. Storage of breast milk

- a. Refrigerate expressed milk immediately. If you express more than once a day, you can add the amounts so long as you put it back in the fridge each time.
- b. Store milk in the fridge for up to 24 hours, or in the freezer for 3 months,
- c. You will notice that stored milk separates out into layers. Shake to remix thoroughly before feeding.
- d. When freezing milk, allow space at the top of the bottle for expansion. It is easiest to freeze milk in individual feed quantities of 60-125ml (2-4ozs)
- e. Thaw frozen milk in the fridge or by standing the bottle in lukewarm water until it has reached body temperature.
- f. If you want to take expressed milk with you when you are out, make sure it stays cold or frozen until needed.

| Method of storing milk       | Use within | REMEMBER  |
|------------------------------|------------|---|
| Refrigerated only            | 24 hours   | <ul style="list-style-type: none"> <li>• Never defrost or heat milk in a microwave</li> <li>• Never refreeze breast milk</li> <li>• Once defrosted, the milk should be kept in the fridge and used within 12 hours</li> <li>• Do not boil breast milk</li> <li>• Throw away any leftover milk after a feed</li> </ul> |
| Refrigerated before freezing | 24 hours   |   |
| Once frozen                  | 3 months   |   |
| Defrosted before feeding     | 12 hours   |   |
| Refreezing                   | NEVER      |   |

#### **K. USING BREAST MILK STORAGE SYSTEM AS A SIPPY CUP**

1. If milk has been chilled, remove lid and place in warm water to heat to temperature required. Always check food temperature before feeding.
2. Hold cup up to baby's mouth and gently tip to allow baby to sip from cup.

**For your child's safety and health  
WARNING!**

**Always use this product with adult supervision. Accidents have occurred when babies have been left alone with drinking equipment due to the baby falling or if the product has disassembled.**

Always check food temperature before feeding. Heating in a microwave oven may produce localised high temperatures. Take extra care when microwave heating. Always stir heated food to ensure even heat distribution and test the temperature before serving.

#### **Cleaning and sterilising**

**Before first use, clean the product**

**Clean before each use.**

**Always wash the product immediately after use**

**Wash in warm soapy water and rinse in clean water or use a dishwasher. Do not clean with solvents or harsh chemicals.**

**Do not allow product to come into contact with oil based foods (e.g. oil/tomato based sauces) as staining will occur**

**Suitable for all forms of sterilisation**

#### **Storage and Use**

**Inspect before each use. Throw away at the first signs of damage or weakness.**

**Do not store or leave in direct sunlight or near a source of heat.  
Do not store with or allow to come into contact with solvents or harsh chemicals  
Not suitable for use in a conventional oven or a combination microwave oven with the  
grill feature switched on.  
Keep all components not in use out of reach of children**

Mayborn Baby and Child  
Northumberland Business Park West  
Cramlington  
Northumberland  
NE23 7RH, UK  
**(or distributor contact details)**

Mayborn Baby and Child is a trading style of Jackel International Limited.  
Copyright© Jackel International Limited 2009

Jackel International Limited retains all intellectual property rights to this product, packaging, and printed or electronic media, including all two and three dimensional designs, graphic images, logos and text.  
Any manufacture, publication, copying, hiring, lending and reproduction of any materials is strictly prohibited and constitutes infringement. TOMMEE TIPPEE and CLOSER TO NATURE are registered trademarks of Jackel International Limited,

Leaflet no. 42341481